Investigating the solvated electron in transient solvent structures utilizing pump-repump-probe spectroscopy

A. Thaller, R. Laenen and A. Laubereau

Physik-Department E11 Technical University of Munich

15.03.2005

- first observation published 140 years ago
- still of high interest:
 - important role in biochemical processes
 - tiny probe for H-bonded systems: Water, methanol, etc.
- generation: electrochemically, radiolysis, photolysis
- generation dynamics occurs on ps-timescale
- first pump-probe-measurements of equilibrated e_{solv} by Barbara et al.
- still open questions regarding precursing states
- $\rightarrow\,$ spectroscopy of transient species $\underline{in\ methanol}$



The solvated electron Properties of e_{solv} :CH₃OH























Generation Process





Generation Process



Generation Process

Excitation of the Equilibrated Groundstate Excitation of the Hot Groundstate Excitation of the Early Bound States Microscopic Model

E11



Generation Process

Excitation of the Equilibrated Groundstate Excitation of the Hot Groundstate Excitation of the Early Bound States Microscopic Model

E11



Generation Process



Generation Process



































